### **Congratulations on your Veet InfiniSilk™ Pro purchase.**

As Veet InfiniSilk™ Pro is a powerful device, it should be used with special attention to safety. Please read the entire user manual before using the product, and adhere to the instructions regarding additional conditions that may warrant safety precautions and warnings.

- Only use this product on lighter skin tones. Do not use Veet InfiniSilk™ Pro on naturally dark skin.
- Do not use on tanned skin or immediately after sun exposure.
- Do not use on permanent make-up or tattooed skin areas.
- Keep away from your eyes!

Note! Veet InfiniSilk™ Pro may not be effective on naturally white and light hair.

## **Quick Start**

### **Plan your Veet InfiniSilk™ Pro hair removal treatments for best results.**

Treatments 1-4 – plan two weeks apart

Treatments 5-7 – plan four weeks apart

Treatments 8 + – treat as needed, until desired results are achieved.

### **Set up Veet InfiniSilk™ Pro!**

1. Remove Veet InfiniSilk™ Pro Pulse Light Technology Device and other components from the box.
2. Verify that the LAMP CARTRIDGE for legs and body is inserted properly at the tip of the device. If you would like to use the smaller, precision lamp cartridge for use on upper lip, cheek, or chin, follow the replacement instructions.
3. Plug the ADAPTOR cord into the Veet InfiniSilk™ Pro HAIR REMOVAL DEVICE socket.
4. Plug the ADAPTOR into an electrical outlet. The STAND BY INDICATOR LIGHT will turn on.

Your Veet InfiniSilk™ Pro is now ready to start.

**Treat with Veet InfiniSilk™ Pro safely and properly!**

5. Choose areas to treat that have not had recent exposure to the sun.
6. The skin of the treatment areas should be shaved, clean, dry and free of any powders, antiperspirants or deodorants. Do not wax, pluck or tweeze the hair out before using the Veet InfiniSilk™ Pro device. Shaving is an important step to get the results you want!
7. Press the ON/STBY SWITCH. The POWER ON INDICATOR LIGHT (orange) will turn on and a fan sound (similar to the sound of a hairdryer) will start.
8. Approximately 3.5 seconds after pressing the ON/STBY SWITCH, the two READY INDICATOR LIGHTS (green) will turn on. The device is then ready for you to trigger the first pulse at the lowest energy level.
9. Before each treatment, choose the desired energy level.  
Veet InfiniSilk™ Pro offers 5 levels of energy from 1 (the lowest) to 5 (the highest).  
To choose the right energy level for your first treatment, test each treatment area separately:

**Skin Test**

- Place the TREATMENT SURFACE on the skin and press the trigger to deliver one pulse.
- If you did not experience abnormal discomfort press the + button once and deliver one pulse on a different spot at energy level 2.
- Continue testing at increased energy levels until you have reached the highest level that you feel comfortable with. If within one hour you experience no adverse effects, you can start your full treatment at that energy level.
- Repeat this test for each body part you intend to treat.

10. Set the energy level. If you want to change the energy level setting, press the ENERGY LEVEL SETTING BUTTONS using the “-“ or “+“ to respectively decrease or increase the energy level, until the desired energy level is set and represented by the energy level indicator lights (For details see “Energy Level” section inside this user manual).

11. Apply the TREATMENT SURFACE to the skin, making sure the skin is spread evenly and smoothly. As soon as the TREATMENT SURFACE is in full contact with the skin, the READY INDICATOR LIGHTS will start blinking.
12. Press the PULSE BUTTON. The device will first determine the colour of your skin. If the colour of the skin is light enough for safe application, the device will flash a pulse of light onto your skin, and the READY INDICATOR LIGHTS will turn off. You will see a bright flash of light and simultaneously hear a subtle pop sound, which is a normal noise for the device. You will feel a mild sensation of warmth and tingling. Veet InfiniSilk™ Pro will immediately recharge for the next pulse. After 3.5 seconds the READY INDICATOR LIGHTS will turn on again.  
NOTE : As the PRECISION LAMP CARTRIDGE does not include a SKIN COLOUR SENSOR, you must first confirm that the area to be treated passes the Skin Colour Sensor test by pulsing it with the larger cartridge.
13. Remove the TREATMENT SURFACE from the treated area of skin.
14. Move the TREATMENT SURFACE to another area of skin. Use the pressure marks the device just made on your skin to guide you for proper positioning of the next pulse, avoiding both gaps and overlaps between pulses.

Try to avoid overlapping pulses!

Warning: Do not treat the same area of skin more than once per hair removal session! Treating the same area of skin more than once per session increases the likelihood of adverse effects.

15. Repeat the process until the entire treatment area is covered, without overlapping pulses.

If your skin blisters or burns, STOP USE IMMEDIATELY!

Copyright © 2014 RB. All rights reserved.  
Print date: August 2014

UK - PO Box 4044, Slough SL1 0NS  
RO: Citywest Business Campus, Dublin 24  
For questions or concerns, call :  
UK – 0845 769 7079 RO: 01 661 7318  
www.veet.co.uk

VEET, VEET logo and Veet INFINI'SILK are trademarks of the Reckitt Benckiser Group of Companies.  
Silk'n and Silk'n logo are trademarks of Home Skinovations Ltd, used under license.

## Table of Contents

<b>1. Before you start</b> .....	5
1.1 Veet Infini'Silk™ Pro device description .....	5
1.2. Package contents .....	6
1.3. The LAMP CARTRIDGE .....	6
<b>2. Veet Infini'Silk™ Pro intended use</b> .....	6
<b>3. Safety with Veet Infini'Silk™ Pro</b> .....	7
<b>4. Contraindications</b> .....	9
<b>5. Avoid using Veet Infini'Silk™ Pro if any of the following applies:</b> .....	12
<b>6. Do not use Veet Infini'Silk™ Pro if:</b> .....	12
<b>7. Possible side effects using Veet Infini'Silk™ Pro</b> .....	13
<b>8. Permanent Hair Reduction with Veet Infini'Silk™ Pro</b> .....	14
<b>9. Get started!</b> .....	15
9.1. What to expect when treating with Veet Infini'Silk™ Pro? .....	19
9.2. After treating with Veet Infini'Silk™ Pro .....	19
<b>10. Maintenance of Veet Infini'Silk™ Pro</b> .....	20
10.1. Cleaning Veet Infini'Silk™ Pro .....	20
10.2. Replacing the LAMP CARTRIDGE .....	20
<b>11. Troubleshooting</b> .....	22
11.1. "My Veet Infini'Silk™ Pro does not start.".....	22
11.2. "A light pulse is not emitted when I press the PULSE BUTTON".....	22
<b>12. Consumer Relations</b> .....	23
<b>13. Specifications</b> .....	23
<b>14. Labelling</b> .....	24

Before using Veet InfiniSilk™ Pro for the first time, please read this User Manual in its entirety. Pay particular attention to sections on device use procedures, device operation, and after-use procedures. We recommend you re-familiarise yourself with this User Manual before each use of Veet InfiniSilk™ Pro. Veet InfiniSilk™ Pro is a powerful electrical device. As such, it should be used with special attention to safety. Please read all warnings and safety precautions before use, and strictly follow them when using Veet InfiniSilk™ Pro.

## 1. Before you start

Veet InfiniSilk™ Pro is a light-based device for permanent hair reduction\* results, designed for home-use.

### 1.1 Veet InfiniSilk™ Pro device description

Your Veet InfiniSilk™ Pro consists of a HAIR REMOVAL DEVICE with a CONTROL PANEL and a LAMP CARTRIDGE. On the Veet InfiniSilk™ Pro HAIR REMOVAL DEVICE HANDLE you can find the PULSE BUTTON. The LAMP CARTRIDGE is located at the tip of the HAIR REMOVAL DEVICE.

- ① Hair Removal device
- ② Control panel
- ③ Pulse button
- ④ Lamp cartridge
- ⑤ Skin Colour Sensor
- ⑥ Treatment Surface
- ⑦ Control Panel
- ⑧ ON/STBY switch
- ⑨ Stand by indicator light
- ⑩ Power on indicator light
- ⑪ Energy level setting buttons
- ⑫ Energy level indicator lights
- ⑬ Replace cartridge indicator light
- ⑭ System warning indicator light
- ⑮ Skin tone warning indicator light
- ⑯ Ready indicator light

\* average hair reduction 54% at 3 months, and 50 % at 18 months post last treatment (study on 62 people)

## 1.2. Package contents

Upon opening the Veet Infini'Silk™ Pro package, you will find the following parts:

- Veet Infini'Silk™ Pro Hair Removal Device with lamp cartridge for legs and body
- A smaller, precision lamp cartridge for usage on upper lip, cheek and chin
- A wall mount adaptor
- This User Manual

## 1.3. The LAMP CARTRIDGE

The Veet Infini'Silk™ Pro LAMP CARTRIDGE can fire a limited and pre-set number of light pulses. The LAMP CARTRIDGE with this device delivers up to 120,000 pulses. In addition to the regular LAMP CARTRIDGE, a Precision lamp cartridge for small and hard to reach areas is also included.

## 2. Veet Infini'Silk™ Pro intended use

**Veet Infini'Silk™ Pro is intended to affect long term, or permanent hair reduction. Permanent hair reduction is defined as a long-term stable reduction in the number of hairs re-growing after a treatment regime.**

Ideal body areas for Veet Infini'Silk™ Pro use include the underarms, bikini line, stomach, arms, legs, back, shoulders, chest, upper lip, cheek and chin.

### 3. Safety with Veet Infini'Silk™ Pro

#### **Important Safety Information – Read Before Use!**

Read this User Manual carefully before you use the appliance and keep it for future reference.

- Keep out of reach of children
- This appliance is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.

**DANGER:** As with most electrical appliances, electrical parts are electrically live even when unplugged.

To reduce the risk of electric shock:

- Do not use near or in bathtubs, showers, sinks or other vessels containing water.
- Do not use with wet hands.
- Do not place or store the appliance where it can fall into a tub or sink.
- Do not place in or drop into water or other liquid.
- If the appliance falls into water while plugged in, unplug it immediately. Do not reach into the water.

**WARNING:** Should this product malfunction, do not attempt to repair it.

This appliance has no user-serviceable parts or components. If the main cord is damaged, you must contact RB Consumer Relations to avoid hazard.

To reduce the risk of burns, electrocution, fire or injury to persons:

- Check if the mains voltage indicated on the appliance corresponds to your local mains voltage before you connect the appliance. The use of adapters, multiple sockets and/or extension cables is not permitted.
- The appliance should never be left unattended when plugged in.
- Use this appliance only for its intended use as described in this leaflet. Only use the power lead included in this pack.
- Never operate this appliance if it has a damaged power lead, if it is not working properly, or if it has been dropped, damaged or dropped in water.
- Keep the power lead away from heated surfaces. Do not wrap the power lead around the appliance.
- Always unplug the appliance after use. Wait until the appliance has cooled down before you store it.

Opening the Veet Infini'Silk™ Pro may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury. It may also damage your device. Please contact RB Consumer Relations if you have a broken or damaged device.

With Veet Infini'Silk™ Pro safety comes first.

- **Veet Infini'Silk™ Pro - Sensor for the skin**

Light based hair removal is not suitable for naturally darker skin tones or tanned skin. Veet Infini'Silk™ Pro comes with a built-in SKIN COLOUR SENSOR that is designed to measure the complexion of the applied surface and enable application only on suitable skin complexions. This safety feature will not let you treat where your skin is too dark or too tanned.

Furthermore, the 1.3X3cm<sup>2</sup> TREATMENT SURFACE through which pulses of light are delivered is recessed inside the LAMP CARTRIDGE. This enables Veet Infini'Silk™ Pro to protect your skin by avoiding direct contact between the TREATMENT SURFACE and the skin.

- **Veet Infini'Silk™ Pro - no need for eyewear**

Veet Infini'Silk™ Pro has a built-in safety feature for eye protection. It has been designed so that a light pulse can not be emitted when the TREATMENT SURFACE is facing open air. The safety switch enables treatment only when the TREATMENT SURFACE is in contact with the skin. It is recommended not to look directly at the flash emitted by the lamp during the treatment.

## 4. Contraindications

### Important safety information – read before use!

Veet Infini'Silk™ Pro is not designed for everyone. Please read and consider the information in the following section before use.

Hair Colour	Skintone					
	White	Beige	Light Brown	Mid Brown	Dark Brown	Black
Black	✓	✓	✓	✓	✗	✗
Dark Brown	✓	✓	✓	✓	✗	✗
Brown	✓	✓	✓	✓	✗	✗
Dark Blond/Light Brown	✓	✓	✓	✓	✗	✗
White/Grey/Red/Light blond	✗	✗	✗	✗	✗	✗

✓ This product is suitable for you to use

✗ This product is neither safe nor suitable for your skintone and hair type

Warning	Possible results	Safety tip
<p><b>Do not use Veet InfiniSilk™ Pro on naturally dark skin complexion.</b></p> <p style="text-align: right;"><b>Not Safe!</b></p> 	<p>Treating dark skin with Veet InfiniSilk™ Pro can result in adverse effects such as redness and skin discomfort</p>	<p>A SKIN COLOUR SENSOR is embedded in Veet InfiniSilk™ Pro to measure the treated skin complexion at the beginning of each session and occasionally during the session. The SKIN COLOUR SENSOR ensures that pulses will only be emitted on suitable skin tones.</p>
<p><b>Do not use Veet InfiniSilk™ Pro on tanned skin or after recent sun exposure (including artificial tanning machines). Do not expose treated areas to the sun.</b></p>	<p>Treating with Veet InfiniSilk™ Pro before or after any sun exposure can result in adverse effects such as redness and skin discomfort. Tanned skin, particularly following sun exposure, contains large quantities of the pigment Melanin. This applies to all skin types and complexions, including those which don't seem to tan quickly. The presence of large quantities of Melanin exposes the skin to higher risk when using any method of light based hair removal.</p>	<p>Avoid exposure to the sun for 4 weeks before your Veet InfiniSilk™ Pro treatment!</p> <p>Sun exposure includes unprotected exposure to direct sunlight of over 15 minutes constantly, or unprotected exposure to diffused sunlight of over 1 hour constantly. To protect recently treated skin when exposed to sunlight, be sure to thoroughly apply sunscreen SPF 30 or higher, for 2 weeks after each hair removal session.</p>

Warning	Possible results	Safety tip
Never use Veet InfiniSilk™ Pro around or near the eyes.		Veet InfiniSilk™ Pro has a built-in safety feature for eye protection. It has been designed so that a light pulse can not be emitted when the TREATMENT SURFACE is facing open air. The safety switch enables treatment only when the TREATMENT SURFACE is in contact with the skin.
Do not treat on tattoos or permanent makeup, dark brown or black spots, (such as large freckles, birth marks, moles or warts), nipples, genitals or lips.	Treating any area with the Veet InfiniSilk™ Pro that is dark in colour or has more pigment can result in adverse effects such as burns, blisters, and skin colour changes (hyper- or hypopigmentation).	Cover area with material that will not absorb the light such as a white cloth or white medical tape.
Do not treat on active eczema, psoriasis, lesions, open wounds or infections (cold sores), abnormal skin conditions caused by systemic or metabolic diseases (diabetes for example)	Light based treatments can cause adverse reactions on already sensitive areas.	Wait for the affected area to heal before using Veet InfiniSilk™ Pro.
A treatment area that has a history of herpes outbreaks.	Light based treatments can cause adverse reactions on already sensitive areas.	Consult your physician to receive preventative treatment before using Veet InfiniSilk™ Pro.
Do not treat with Veet InfiniSilk™ Pro if you are pregnant or nursing.	Hormonal changes during pregnancy could increase skin sensitivity and the risk of adverse effects.	
Do not treat on areas where you may want hair to grow.	Hair reduction results are permanent.	

## **5. Avoid using Veet Infini'Silk™ Pro if any of the following applies:**

- If you have a history of keloidal scar formation.
- If you suffer from epilepsy.
- If you have a disease related to photosensitivity, such as porphyria, polymorphic light eruption, solar urticaria, lupus, etc.
- If you have a history of skin cancer or areas of potential skin malignancies.
- If you have received radiation therapy or chemotherapy treatments within the past 3 months.
- If you have any other condition which in your physician's opinion would make it unsafe for you to be treated.
- If your skin was treated with a physician-prescribed drug within the past 6 months, please consult with your physician.

If any of these apply to you, please consult with your physician prior to using the Veet Infini'Silk™ Pro. If you have a doubt regarding the use of this device, please contact your physician.

## **6. Do not use Veet Infini'Silk™ Pro if:**

- the LAMP CARTRIDGE is damaged, or if its TREATMENT SURFACE is cracked, chipped or missing.
- the attachments or accessories are not recommended by RB.
- it has a damaged cord or plug. Keep the power cord away from heated surfaces.
- you see or smell smoke when it is in use.
- it is not working properly or if it appears damaged.
- its SKIN COLOUR SENSOR is cracked, or broken.
- the outer shell is cracked or is coming apart.

If any of these have happened to your device, stop use immediately and contact RB Consumer Relations.

## 7. Possible side effects using Veet InfiniSilk™ Pro

When used according to the instructions, side effects and complications associated with use of Veet InfiniSilk™ Pro are uncommon. However every cosmetic procedure, including those designed for home use, involves some degree of risk. Therefore it is important that you understand and accept the risks and complications that can occur with pulsed light hair removal systems designed for home use.

Adverse reaction	Degree of risk
<p><b>Minor skin discomfort</b> – Although home pulsed light hair removal is generally very well-tolerated, most users do feel some mild discomfort during use, usually described as being a mild stinging sensation on the treated skin areas. The stinging sensation usually lasts during the time of the application itself or for a few minutes thereafter. Anything beyond this minor discomfort is abnormal and means that either you should not continue to use Veet InfiniSilk™ Pro because you are unable to tolerate the hair removal application, or that the energy level setting is too high.</p>	Minor
<p><b>Skin redness</b> – Your skin may become red right after using Veet InfiniSilk™ Pro or within 24 hours of using Veet InfiniSilk™ Pro. Redness generally clears up within 24 hours. See your doctor if redness does not go away within 2 to 3 days.</p>	Minor
<p><b>Increased sensitivity to skin</b> – The skin of the treated area is more sensitive so you may encounter dryness or flaking of the skin.</p>	Minor
<p><b>Pigment changes</b> – Veet InfiniSilk™ Pro targets the hair shaft, in particular the pigmented cells in the hair follicle and the hair follicle itself. Nevertheless there is risk of temporary hyperpigmentation (increased pigment or brown discoloration) or hypopigmentation (whitening) to the surrounding skin. This risk of changes in skin pigmentation is higher for people with darker skin tones. Usually discoloration or changes to skin pigment are temporary and permanent hyperpigmentation or hypopigmentation rarely occur.</p>	Rare
<p><b>Excessive redness and swelling</b> – In rare cases treated skin may become very red and swollen. This is more common in sensitive areas of the body. The redness and swelling should subside within 2 to 7 days and should be treated with frequent applications of ice. Gentle cleansing is OK, but you should avoid exposure to the sun.</p>	Rare

## 8. Permanent hair reduction with Veet Infini'Silk™ Pro

A typical full hair growth cycle may take 18-24 months. During this time multiple Veet Infini'Silk™ Pro sessions may be required in order to achieve permanent hair reduction. The efficiency of hair removal varies from person to person according to body area, hair colour, and how Veet Infini'Silk™ Pro is used.

### Typical Veet Infini'Silk™ Pro hair removal plan during a full hair growth cycle:

Treatments 1-4 – plan two weeks apart

Treatments 5-7 – plan four weeks apart

Treatments 8 + – treat as needed, until desired results are achieved.

### Expected treatment outcome for hair removal

Fitzpatrick skin type	Hair Colour	Light Fluence [J/cm <sup>2</sup> ]	Anatomic area of the body <sup>2</sup>	Average Number of Treatments <sup>1</sup>	Effectiveness (% of hair removal/regrowth) <sup>1</sup>
I-IV	Light brown to brown	3-4.5	Lower or upper legs	Session 10-12, approx. four weeks apart.	60% Hair reduction
			Arms		50% Hair reduction
			Bikini Line		60% Hair reduction
			Axilla		55% Hair reduction
I-IV	Brown to black	3-4.5	Lower or upper legs	Session 8-10, approx. four weeks apart.	70% Hair reduction
			Arms		60% Hair reduction
			Bikini Line		70% Hair reduction
			Axilla		65% Hair reduction

- Individual response depends on hair type as well as biological factors that may affect hair growth patterns. Some users may respond quicker or slower than the average number of treatments.
- One cannot expect permanent hair reduction in a single or even double treatment sessions. The duration of the resting period for hair follicles also depends on the body area.

## 9. Get started!

1. Remove Veet Infini'Silk™ Pro HAIR REMOVAL device and other components from box.
2. Verify that the LAMP CARTRIDGE for legs and body is inserted properly at the tip of the device. If you would like to use the smaller, precision lamp cartridge for use on upper lip, cheek, or chin, follow the replacement instructions.
3. Plug the ADAPTOR cord into the Veet Infini'Silk™ Pro HAIR REMOVAL DEVICE socket.
4. Plug the ADAPTOR into an electrical outlet. The STAND BY INDICATOR LIGHT will turn on.

Your Veet Infini'Silk™ Pro is now ready to start.

**Treat with Veet Infini'Silk™ Pro safely and properly!**  
**Choose areas to treat that have not had recent exposure to the sun.**

5. The skin of the treatment areas should be shaved, clean, dry and free of any powders, antiperspirants or deodorants. Do not wax, pluck or tweeze the hair out before using the Veet Infini'Silk™ Pro device. Shaving is an important step to get the results you want!
6. Press the ON/STBY SWITCH. The POWER ON INDICATOR LIGHT (orange) will turn on and a fan sound (similar to the sound of a hairdryer) will start.
7. Approximately 3.5 seconds after pressing the ON/ STBY SWITCH, the two READY INDICATOR LIGHTS (green) will turn on. The device is then ready for you to trigger the first pulse at the lowest energy level.
8. Before each treatment, choose the energy level. Veet Infini'Silk™ Pro offers 5 levels of energy from 1 (the lowest) to 5 (the highest). To choose the right energy level for your first treatment, test each treatment area separately:

### Skin Test

- Place the TREATMENT SURFACE on the skin and press the trigger to deliver one pulse. If you did not experience abnormal discomfort press the + button once and deliver one pulse on a different spot at energy level 2.
- Continue testing at increased energy levels until you have reached the highest level that you feel comfortable with. If within one hour you experience no adverse effects, you can start your full treatment at that energy level.
- Repeat this test for each body part you intend to treat.

9. Set the energy level. If you want to change the energy level setting, press the ENERGY LEVEL SETTING BUTTONS using the “-“ or “+“ to respectively decrease or increase the energy level, until the desired energy level is set and represented by the energy level indicator lights (For details see “Energy Level” section inside this user manual).
10. Apply the TREATMENT SURFACE to the skin, making sure the skin is spread evenly and smoothly. As soon as the TREATMENT SURFACE is in full contact with the skin, the READY INDICATOR LIGHTS will start blinking.
11. Press the PULSE BUTTON. The device will first determine the colour of your skin. If the colour of the skin is light enough for safe application, the device will flash a pulse of light onto your skin, and the READY INDICATOR LIGHTS will turn off. You will see a bright flash of light and simultaneously hear a subtle pop sound, which is a normal noise for the device. You will feel a mild sensation of warmth and tingling. Veet InfiniSilk™ Pro will immediately recharge for the next pulse. After 3.5 seconds the READY INDICATOR LIGHTS will turn on again.  
NOTE : As the PRECISION LAMP CARTRIDGE does not include a SKIN COLOUR SENSOR, you must first confirm that the area to be treated passes the Skin Colour Sensor test by pulsing it with the larger cartridge.
12. Remove the TREATMENT SURFACE from the treated area of skin.
13. Move the TREATMENT SURFACE to another area of skin. Use the pressure marks the device just made on your skin to guide you for proper positioning of the next pulse, avoiding both gaps and overlaps between pulses.

Veet InfiniSilk™ Pro pulses should be administered in rows, starting at one end of each row and progressing sequentially towards the other end. This technique allows better control of skin coverage, and helps you avoid treating the same area more than once or overlapping skin areas.

When applied to the skin, the Veet InfiniSilk™ Pro TREATMENT SURFACE is designed to create temporary pressure marks on the treated area. These visible marks can be used for proper positioning of the next pulse.

Try to avoid overlapping pulses!

**Warning:** Do not treat the same area of skin more than once per hair removal session! Treating the same area of skin more than once per session increases the likelihood of adverse effects.

14. Repeat the process until the entire treatment area is covered, without overlapping pulses.

If your skin blisters or burns, STOP USE IMMEDIATELY!

### **Energy level**

Energy level determines the intensity of the Veet InfiniSilk™ Pro light pulse delivered to your skin, from the lowest level (-) to the highest level (+). As energy level increases, so does hair removal results as well as the risk of possible side effects and complications.

Veet InfiniSilk™ Pro offers 5 levels of energy from 1 (the lowest) to 5 (the highest), represented by 3 energy level indicator lights:

Level 1 (lowest) - the left indicator light is steadily lit

Level 2 - the left indicator light is steadily lit and the central indicator light is blinking.

Level 3 - both left and central indicator lights are steadily lit

Level 4 - the left and central indicator lights are steadily lit and the right indicator light is blinking.

Level 5 (highest) - all 3 indicator lights are steadily lit.

Whenever Veet InfiniSilk™ Pro is turned on its energy level will automatically be reset to the lowest energy level. Only one ENERGY LEVEL INDICATOR LIGHT will be on.

To set the energy level, press the ENERGY LEVEL SETTING BUTTONS using “-“ or “+” to respectively decrease or increase the energy level.

### **The skin colour sensor**

Light-based hair removal on darker skin tones can result in adverse effects such as burns, blisters, and skin colour changes (hyper- or hypo-pigmentation).

To prevent such misuse, a SKIN COLOUR SENSOR in Veet InfiniSilk™ Pro measures the treated skin complexion at the beginning of each session and occasionally during the session. If the SKIN COLOUR SENSOR detects a skin tone that is too dark for Veet InfiniSilk™ Pro application, the device will automatically stop emitting pulses. If you see no light pulse and the SKIN TONE WARNING INDICATOR LIGHT blinks, this is an indication that your skin tone, measured by the SKIN COLOUR SENSOR, is too dark for safe application. Try using the device on a different body part or contact RB Consumer Relations.

NOTE : As the PRECISION LAMP CARTRIDGE does not include a SKIN COLOUR SENSOR, you must first confirm that the area to be treated passes the Skin Colour Sensor test by pulsing it with the larger cartridge.

### 9.1. What to expect when treating with Veet Infini'Silk™ Pro?

For many people, using Veet Infini'Silk™ Pro may be their first experience with a light-based device designed for home use. Veet Infini'Silk™ Pro is simple to use, and hair removal sessions go by quickly. During a Veet Infini'Silk™ Pro session it is normal to experience and feel:

- **A fan noise** – The cooling fan in Veet Infini'Silk™ Pro makes noise similar to a hairdryer. This is normal.
- **A pop sound with each pulse** – When a pulse of light is activated, it is normal to hear a subtle pop sound simultaneously with the flash of light.
- **Moderate pressure of the treatment surface** – This is necessary and helpful for placement of adjacent pulses of light, and is part of the safety feature of Veet Infini'Silk™ Pro.
- **A sensation of warmth and tingling** – During each pulse of light it is normal to feel a mild sensation of warmth and tingling from the light energy.
- **Some mild red or pink colour** – During and just after your Veet Infini'Silk™ Pro session it is not uncommon to see some very mild, pink-like colour of the skin. This is usually most noticeable around the hairs themselves.

However if you see full redness of the skin, blistering or burns stop use of Veet Infini'Silk™ Pro immediately.

### 9.2. After treating with Veet Infini'Silk™ Pro

- When Veet Infini'Silk™ Pro session has been completed turn Veet Infini'Silk™ Pro off by pressing the ON/STBY SWITCH. All the lit indicators will turn off and the STBY INDICATOR LIGHT will turn on (be sure to remember the last energy level setting you used, as it will not be restored when turning Veet Infini'Silk™ Pro on again.)
- Unplug the ADAPTOR cord from the electrical outlet.
- After each hair removal session it is recommended that you clean your Veet Infini'Silk™ Pro device, especially the TREATMENT SURFACE (See: "Cleaning Veet Infini'Silk™ Pro").
- After cleaning, it is recommended to store your Veet Infini'Silk™ Pro device in its original box, and keep it away from water.

**Skin care following hair removal session**

Do not expose treated areas of skin to the sun. Be sure to carefully protect the treated skin with sunscreen, throughout the hair removal period and for at least 2 weeks following the last Veet Infini'Silk™ Pro session.

**10. Maintenance of Veet Infini'Silk™ Pro****10.1. Cleaning Veet Infini'Silk™ Pro**

After each hair removal session, it is recommended to clean your Veet Infini'Silk™ Pro device, and especially the TREATMENT SURFACE.

- Unplug Veet Infini'Silk™ Pro before cleaning.
- Use a dry, clean cloth and a specially formulated cleaner for electronic equipment to gently wipe Veet Infini'Silk™ Pro surface, and especially the TREATMENT SURFACE. Never immerse Veet Infini'Silk™ Pro or any of its parts in water!

**10.2. Replacing the LAMP CARTRIDGE**

The Veet Infini'Silk™ Pro LAMP CARTRIDGE can fire a pre-set and limited number of light pulses. Pulse intensity is determined only according to the energy level setting of the device. There is no decrease of energy during the usable lifetime of the LAMP CARTRIDGE. When only 100 light pulses remain in the LAMP CARTRIDGE, the REPLACE CARTRIDGE INDICATOR LIGHT on the CONTROL PANEL will start blinking, indicating that the LAMP CARTRIDGE should soon be replaced. As the number of light pulses in the LAMP CARTRIDGE further decreases, the REPLACE CARTRIDGE INDICATOR LIGHT blinking will grow faster. Once all the pulses in a LAMP CARTRIDGE have been used, the REPLACE CARTRIDGE INDICATOR LIGHT will remain steadily on, indicating that further pulses cannot be emitted.

In order to continue the hair removal session, it will now be necessary to replace the LAMP CARTRIDGE.

**Note:** A LAMP CARTRIDGE should also be replaced if large spots appear inside it or if its TREATMENT SURFACE is broken.

To replace a LAMP CARTRIDGE:

1. Press the ON/STBY SWITCH to turn Veet InfiniSilk™ Pro to STBY mode.
2. Unplug the power adaptor from the electrical outlet.
3. Grasp the used LAMP CARTRIDGE on both sides, pull it out of the socket.
4. Unwrap a new LAMP CARTRIDGE or the PRECISION LAMP CARTRIDGE.
5. Push the new LAMP CARTRIDGE or the PRECISION LAMP CARTRIDGE gently into place.

Note! After replacing the lamp cartridge with a new one, always set the energy to be one level lower than your normal setting. Raise the energy level back, only if you do not experience abnormal discomfort during or after the treatment.

The PRECISION LAMP CARTRIDGE contains light pulses. With its smaller treatment surface, the PRECISION LAMP CARTRIDGE is designed for smaller body areas and the face. As the PRECISION LAMP CARTRIDGE does not include a SKIN COLOUR SENSOR, you must first confirm that the area to be treated passes the Skin Colour Sensor test by pulsing it with the larger cartridge.

## 11. Troubleshooting

### 11.1. “My Veet Infini’Silk™ Pro does not start.”

Make sure the POWER ADAPTOR is properly connected to the Veet Infini’Silk™ Pro device. Make sure the POWER ADAPTOR is plugged into an electrical outlet on the wall.

### 11.2. “A light pulse is not emitted when I press the PULSE BUTTON”

Make sure that the READY INDICATOR LIGHT is blinking to ensure you have good contact with the skin and that the TREATMENT SURFACE is evenly and firmly pressed to the skin. For your safety, the PULSE BUTTON will activate a pulse only if the TREATMENT SURFACE is firmly pressed against the skin.

- Check the SKIN TONE WARNING INDICATOR LIGHT. If it blinks, this is an indication that your skin tone, measured by the SKIN COLOUR SENSOR, is too dark for safe application. Try using the device on a different body part or contact RB Consumer Relations.
- Check the REPLACE CARTRIDGE INDICATOR LIGHT. If it is on, disconnect Veet Infini’Silk™ Pro and replace the LAMP CARTRIDGE.
- Make sure that the READY INDICATOR is on.
  - If within 10 seconds the READY INDICATOR remains off turn Veet Infini’Silk™ Pro off and back on by pressing the ON/STBY SWITCH twice.
- Check the SYSTEM WARNING INDICATOR LIGHT.
  - If the light is on turn Veet Infini’Silk™ Pro off and back on by pressing the ON/STBY SWITCH twice.

If the problem persists, contact RB Consumer Relations UK – 0845 769 7079  
ROI: 01 661 7318 or visit [www.veet.co.uk](http://www.veet.co.uk)

Do not attempt to open or repair your Veet Infini’Silk™ Pro device.

## 12. Consumer Relations

For more information about Veet InfiniSilk™ Pro or Veet products, contact the RB Consumer Relations Centre - UK: 0845 769 7079 ROI: 01 661 7318 or visit [www.veet.co.uk](http://www.veet.co.uk)  
If your device is broken or damaged, or for any other user assistance, contact the RB Consumer Relations Centre - UK: 0845 769 7079 ROI: 01 661 7318 or visit [www.veet.co.uk](http://www.veet.co.uk)

## 13. Specifications

Spot Size	1.3cm x 3cm [3.9cm <sup>2</sup> ]	
Speed	1 pulse every 3.5 second: 1.1 cm <sup>2</sup> /sec	
Technology	Home Pulsed Light	
Max Energy Level	Max 5J/cm <sup>2</sup>	
Wavelength	475-1200nm	
Charging Time/Power Source	Continuous operation	
Electrical Requirements	100-240VAC, 0.9A	
Operation and Safety	Skin Colour Sensor ensures use only on suitable skin complexions.	
Package Size	22cm x 22cm x 14cm	
System weight	0.38kg	
Temperature	Operating	10 to 40°C
	Storage	-40 to +70°C
Relative Humidity	Operating	30 to 75%rH
	Storage	10 to 90%rH
Atmospheric pressure	Operating	700 to 1060hPa
	Storage	500 to 1060hPa
Adaptor Rated	Model: KSAPV0361200300DS 100-240V; 50-60Hz; 0.9A	

## 14. Labelling

Degree of protection against electric shock: type BF applied part

Please DO NOT dispose of appliance in normal household waste. Your local appliance retailer may operate a 'take-back' scheme when you are ready to purchase a replacement product, alternatively contact your local government authority for further help and advice on where to take your appliance for recycling.

This product carries the CE mark and conforms to EC Directives.

Follow operating instructions

Degree of protection against ingress water: ordinary

This device is not suitable for use in the presence of flammable anesthetic mixture with air or with Oxygen or Nitrous Oxide.

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Copyright 2014 © RB. Alle Rechte vorbehalten.
Copyright 2014 © RB. Tutti i diritti riservati.